



# T-Cup<sup>®</sup>

## Urine Drug Test



- CLIA Waived, 510(K) Cleared & Forensic Use Only
- Female User Friendly & Customizable
- Test up to 15 Drugs per Cup

### CLIA Waived Drugs

	Amphetamine (AMP)	1000ng/mL 500ng/mL		Tricyclic Antidepressants (TCA)	1000ng/mL
	Barbiturates (BAR)	300ng/mL		Methamphetamine ( MET/mAMP)	1000ng/mL 500ng/mL
	Benzodiazepines (BZO)	300ng/mL		Methylenedioxymethamphetamine (MDMA)	500ng/mL
	Buprenorphine (BUP)	10ng/mL		Morphine (MOP/OPI300)	300ng/mL
	Cocaine (COC)	300ng/mL 150ng/mL		Opiate (OPI)	2000ng/mL
	Cannabinoids (THC)	50ng/mL		Oxycodone (OXY)	100ng/mL
	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300ng/mL		Phencyclidine (PCP)	25ng/mL
	Methadone (MTD)	300ng/mL		Propoxyphene (PPX)	300ng/mL

Built-in Adulteration Parameters: Creatinine (CR), Specific Gravity (SG) & pH. Oxidant (OX) & Nitrite (NI) available upon request.

Distributed By: Blue Grass Drug Screen, Inc.  
www.uscreeningsource.com 1-866-323-7336

# Order Information

Item Code	Description	Complexity
TDOA-154	5-Panel T-Cup; AMP, COC, OPI, PCP, THC	CLIA WAIVED
TDOA-254	5-Panel T-Cup; AMP, COC, mAMP, OPI, THC	CLIA WAIVED
TDOA-254A3	5-Panel T-Cup; AMP, COC, mAMP, OPI, THC with 3 AD	CLIA WAIVED
TDOA-654	5-Panel T-Cup; BZO, COC, mAMP, OPI, THC	CLIA WAIVED
TDOA-264	6-Panel T-Cup; AMP, BZO, COC, mAMP, OPI, THC	CLIA WAIVED
TDOA-264A3	6-Panel T-Cup; AMP, BZO, COC, mAMP, OPI, THC with 3 AD	CLIA WAIVED
TDOA-564	6-Panel T-Cup; BZO, COC, mAMP, OPI, OXY, THC	CLIA WAIVED
TDOA-274	7-Panel T-Cup; AMP, BZO, COC, mAMP, OPI, OXY, THC	CLIA WAIVED
TDOA-3104	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MTD, OPI, PCP, THC	CLIA WAIVED
TDOA-3104A3	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MTD, OPI, PCP, THC with 3 AD	CLIA WAIVED
TDOA-4104	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MTD, OPI, OXY, PCP, THC	CLIA WAIVED
TDOA-7104	10-Panel T-Cup; BAR, BUP, BZO, COC, mAMP, MTD, OPI, OXY, TCA, THC	CLIA WAIVED
TDOA-7104A3	10-Panel T-Cup; BAR, BUP, BZO, COC, mAMP, MTD, OPI, OXY, TCA, THC with 3 AD	CLIA WAIVED
TDOA-8104	10-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, THC	CLIA WAIVED
TDOA-3124	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, MDMA, MTD, OPI, OXY, PCP, TCA, THC	CLIA WAIVED
TDOA-4124	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MTD, OPI, OXY, PCP, PPX, THC	CLIA WAIVED
TDOA-6124	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, THC	CLIA WAIVED
TDOA-6124A3	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, THC with 3 AD	CLIA WAIVED
TDOA-6125	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC	CLIA WAIVED
TDOA-6125A3	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC with 3 AD	CLIA WAIVED
TDOA-7125	12-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC	CLIA WAIVED
TDOA-7125A3	12-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC with 3 AD	CLIA WAIVED
TDOA-2135	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC	CLIA WAIVED
TDOA-2135A3	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC with 3 AD	CLIA WAIVED
TDOA-6135F	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, FTY	Forensic Use Only
TDOA-7135EFT	13-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG, FTY, TRA	Forensic Use Only
TDOA-8135E	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG	Forensic Use Only
TDOA-1144	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC	CLIA WAIVED
TDOA-1144A3	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC with 3 AD	CLIA WAIVED
TDOA-1145A3	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, PPX, TCA, THC with 3 AD	CLIA WAIVED
TDOA-8145EF	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG, FTY	Forensic Use Only
TDOA-9145EFTK	14-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG, FTY, TRA, K2	Forensic Use Only
TDOA-9145A3EFTK	14-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG, FTY, TRA, K2 <small>W/A3</small>	Forensic Use Only